

Patients

app

Therapy Reports



Otto Ebenbichler recommends:

„The list of Therapy Reports can be mapped alphabetically, numerically ascending (A-Z, 1-100) or descending (Z-A, 100-1) or can be mapped individually.

To do this, click on the title, while holding the left mouse button pressed move the title. The system keeps track of any resulting order. At the next system entry it shows automatically the collection sequence favored by you.”

menue icons

- public access
- observer permission
- descending order
- ascending order
- individual order
- close/open menue

command bar

- create new treatment
- print content
- send content as an email
- duplicates treatment
- patient file
- refresh platform
- download manual
- delete all data

date

The date of a new therapy report can be set with a single mouse click. Click "Date", a calendar window opens. Leave the date chooser by making a simple mouse click outside.

treatment

After a new therapy report has been created you can select an already documented treatment. A simple mouse click on a treatment title imports all relevant data. The grey data are modifiable only in the app "Treatments".

user manual

impress

idea & concept:
Mag. Otto Walter Ebenbichler

technical development
Gunharth Randolf


design
DI Michael Ebenbichler

product responsibility
communautic®ltd

Physio Observer® is a member of the
communautic®group

service requests to:
office@physioobserver.com

all rights reserved by O.W.Ebenbichler



observer

With "Observer" you decide whether users with observer permission are „not authorized" (default setting) or „authorized" to see the content data.


modules

This gray bar shows the modules you are allowed to work with. With one simple mouse click on a module title you switch to the desired module. The position at which you leave the module remains stored, that means by going back you will get to the same place where you have left the module.

autor information

Here you receive "log information". It is documented by whom and when the last entry dates (entry links) and by whom and when this data sheet was created (entry right).

sent

„Sent" works automatically. Once you have sent the contents of this page by , the programm generates an entry to whom and when this content was sent.